

Anti-Collagen I/COL1A1 Antibody Picoband® HRP Conjugated

Catalog Number: PB9939-HRP

About Col1a1

Collagen, type I, alpha 1, also known as COL1A1, is a human gene that encodes the major component of type I collagen, the fibrillar collagen found in most connective tissues, including cartilage. This gene is mapped to 17q21.33. And this gene encodes the [pro-alpha1 chains](#) of type I collagen whose triple helix comprises two alpha1 chains and one alpha2 chain. Type I is a fibril-forming collagen found in most connective tissues and is abundant in bone, cornea, dermis and tendon. Mutations in this gene are associated with osteogenesis imperfecta types I-IV, Ehlers-Danlos syndrome type VIIA, Ehlers-Danlos syndrome Classical type, Caffey Disease and idiopathic osteoporosis.

Overview

Product Name	Anti-Collagen I/COL1A1 Antibody Picoband® HRP Conjugated
Reactive Species	Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P11087

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of mouse Collagen I, different from the related human sequence by three amino acids, and identical to the related rat sequence.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.

ELISA, Optimal dilutions should be determined by end users.

59 Publications Citing This Product

1. PubMed ID: 10.1038/s41598-020-67514-4, Investigating the antifibrotic effect of the antiparasitic drug Praziquantel in in vitro and in vivo preclinical models
2. PubMed ID: 30888602, Wei X,Bao Y,Zhan X,Zhang L,Hao G,Zhou J,Chen Q.MiR-200a ameliorates peritoneal fibrosis and functional deterioration in a rat model of peritoneal dialysis.Int Urol Nephrol.2019 May;51(5):889-896.doi:10.1007/s11255-019-02122-4.Epub 2019 Mar 19.PMID:30888602;PMCID:PMC6499761.
3. PubMed ID: 31169440, Li X,Bu X,Yan F,Wang F,Wei D,Yuan J,Zheng W,Su J,Yuan J.Deletion of discoidin domain receptor 2 attenuates renal interstitial fibrosis in a murine unilateral ureteral obstruction model.Ren Fail.2019 Nov;41(1):481-488.doi:10.1080/0886022X.2019.1621759.PMID:31169440;PMCID:PMC6567249.

Visit bosterbio.com/anti-collagen-i-picoband-trade-antibody-pb9939-boster.html to see all 59 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Collagen I/COL1A1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.