

## Anti-Intestinal FABP/FABP2 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9943-Fluoro488

### About FABP2

FABP 2, Fatty acid-binding protein 2, is a protein that in humans is encoded by the FABP2 gene. Using a human cDNA probe, the gene is assigned to chromosome 4 in somatic cell hybrids. FABP 2 gene contains four exons and is an abundant cytosolic protein in small intestine epithelial cells. The FABPs belong to a multigene family with nearly twenty identified members. And FABPs are divided into at least three distinct types, namely the hepatic-, intestinal- and cardiac-type. They form 14-15 kDa proteins and are thought to participate in the uptake, intracellular metabolism and/or transport of long-chain fatty acids. Also, they may be responsible in the modulation of cell growth and proliferation.

### Overview

Product Name	Anti-Intestinal FABP/FABP2 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P12104

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human FABP2/I-FABP, different from the related mouse sequence by seven amino acids, and from the related rat sequence by six amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 10.1038/s41467-021-27421-2, Dynamic regulation of N<sup>6</sup>,2'-O-dimethyladenosine (m<sup>6</sup>Am) in obesity

Visit [bosterbio.com/anti-fabp2-i-fabp-picoband-trade-antibody-pb9943-boster.html](https://bosterbio.com/anti-fabp2-i-fabp-picoband-trade-antibody-pb9943-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Intestinal FABP/FABP2 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.