

## Anti-Intestinal FABP/FABP2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9943-Fluoro647

### About FABP2

FABP 2, Fatty acid-binding protein 2, is a protein that in humans is encoded by the FABP2 gene. Using a human cDNA probe, the gene is assigned to chromosome 4 in somatic cell hybrids. FABP 2 gene contains four exons and is an abundant cytosolic protein in small intestine epithelial cells. The FABPs belong to a multigene family with nearly twenty identified members. And FABPs are divided into at least three distinct types, namely the hepatic-, intestinal- and cardiac-type. They form 14-15 kDa proteins and are thought to participate in the uptake, intracellular metabolism and/or transport of long-chain fatty acids. Also, they may be responsible in the modulation of cell growth and proliferation.

### Overview

Product Name	Anti-Intestinal FABP/FABP2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P12104

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human FABP2/I-FABP, different from the related mouse sequence by seven amino acids, and from the related rat sequence by six amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 10.1038/s41467-021-27421-2, Dynamic regulation of N6,2'-O-dimethyladenosine (m6Am) in obesity

Visit [bosterbio.com/anti-fabp2-i-fabp-picoband-trade-antibody-pb9943-boster.html](https://bosterbio.com/anti-fabp2-i-fabp-picoband-trade-antibody-pb9943-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Intestinal FABP/FABP2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.