

Anti-FGF1 Antibody Picoband® APC Conjugated

Catalog Number: PB9944-APC

About FGF1

Fibroblast growth factor 1 (acidic), also known as FGF1/ECGF/HBGF1, is a human gene which is mapped to 5q31. The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It also acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis.

Overview

| | |
|----------------------|--|
| Product Name | Anti-FGF1 Antibody Picoband® APC Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P05230 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E. coli-derived human FGF1 recombinant protein (Position: F16-D155). Human FGF1 shares 96.4% amino acid (aa) sequence identity with both mouse and rat FGF1. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

1 Publications Citing This Product

1. PubMed ID: 25667662, Zhang L, Wang H, Wang T, Jiang N, Yu P, Chong Y, Fu F. Exp Ther Med. 2015 Mar;9(3):972-976. Epub 2014 Dec 24. Ferulic Acid Ameliorates Nerve Injury Induced By Cerebral Ischemia In Rats.

Visit bosterbio.com/anti-fgf1-picoband-trade-antibody-pb9944-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGF1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.