

## Anti-FGF1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9944-Fluoro550

### About FGF1

Fibroblast growth factor 1 (acidic), also known as FGF1/ECGF/HBGF1, is a human gene which is mapped to 5q31. The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It also acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-FGF1 Antibody Picoband® Fluoro550 Conjugated  |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | P05230   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E. coli-derived human FGF1 recombinant protein (Position: F16-D155). Human FGF1 shares 96.4% amino acid (aa) sequence identity with both mouse and rat FGF1. |
| Cross Reactivity    | No cross-reactivity with other proteins  |
| Isotype             | Rabbit IgG   |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | Fluoro550<br>Excitation Wavelength: 562 nm<br>Emission Wavelength: 576 nm  |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.   |

## 1 Publications Citing This Product

---

1. PubMed ID: 25667662, Zhang L, Wang H, Wang T, Jiang N, Yu P, Chong Y, Fu F. Exp Ther Med. 2015 Mar;9(3):972-976. Epub 2014 Dec 24. Ferulic Acid Ameliorates Nerve Injury Induced By Cerebral Ischemia In Rats.

Visit [bosterbio.com/anti-fgf1-picoband-trade-antibody-pb9944-boster.html](http://bosterbio.com/anti-fgf1-picoband-trade-antibody-pb9944-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGF1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.