

## Anti-GADD45A Antibody Picoband® Biotin Conjugated

Catalog Number: PB9945-Biotin

### About GADD45A

Growth arrest and DNA-damage-inducible protein GADD45 alpha (GADD45A) is a protein that in humans is encoded by the GADD45A gene. This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. It is located on 1p34-p12. Sequence analysis of human and hamster cDNA clones demonstrated that the gene has been highly conserved and encodes a novel 165-amino acid polypeptide. Northern blot analysis detected moderate expression of a 1.4-kb GADD45A transcript in heart, skeletal muscle, and kidney, with little or no expression in brain, placenta, lung, liver, and pancreas. In addition, Gadd45a promotes epigenetic gene activation by repair-mediated DNA demethylation.

### Overview

Product Name	Anti-GADD45A Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P24522

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human GADD45A, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users.

Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.  
ELISA, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GADD45A Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.