

Anti-GADD45A Antibody Picoband® Fluoro488 Conjugated

Catalog Number: PB9945-Fluoro488

About GADD45A

Growth arrest and DNA-damage-inducible protein GADD45 alpha (GADD45A) is a protein that in humans is encoded by the GADD45A gene. This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. It is located on 1p34-p12. Sequence analysis of human and hamster cDNA clones demonstrated that the gene has been highly conserved and encodes a novel 165-amino acid polypeptide. Northern blot analysis detected moderate expression of a 1.4-kb GADD45A transcript in heart, skeletal muscle, and kidney, with little or no expression in brain, placenta, lung, liver, and pancreas. In addition, Gadd45a promotes epigenetic gene activation by repair-mediated DNA demethylation.

Overview

| | |
|----------------------|--|
| Product Name | Anti-GADD45A Antibody Picoband® Fluoro488 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P24522 |

Technical Details

| | |
|------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human GADD45A, identical to the related mouse and rat sequences. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GADD45A Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.