

Anti-GADD45A Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9945-Fluoro550

About GADD45A

Growth arrest and DNA-damage-inducible protein GADD45 alpha (GADD45A) is a protein that in humans is encoded by the GADD45A gene. This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. It is located on 1p34-p12. Sequence analysis of human and hamster cDNA clones demonstrated that the gene has been highly conserved and encodes a novel 165-amino acid polypeptide. Northern blot analysis detected moderate expression of a 1.4-kb GADD45A transcript in heart, skeletal muscle, and kidney, with little or no expression in brain, placenta, lung, liver, and pancreas. In addition, Gadd45a promotes epigenetic gene activation by repair-mediated DNA demethylation.

Overview

Product Name	Anti-GADD45A Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P24522

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human GADD45A, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GADD45A Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.