

Anti-ISG15 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9950-Fluoro647

About ISG15

Interferon-stimulated gene 15 (ISG15) is a 17 kDA secreted protein that in humans is encoded by the ISG15 gene. The protein encoded by this gene is a ubiquitin-like protein that is conjugated to intracellular target proteins upon activation by interferon-alpha and interferon-beta. Several functions have been ascribed to the encoded protein, including chemotactic activity towards neutrophils, direction of ligated target proteins to intermediate filaments, cell-to-cell signaling, and antiviral activity during viral infections. While conjugates of this protein have been found to be noncovalently attached to intermediate filaments, this protein is sometimes secreted.

Overview

| | |
|----------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-ISG15 Antibody Picoband® Fluoro647 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P05161 |

Technical Details

| | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen | E.coli-derived human ISG15 recombinant protein (Position: G2-G157). Human ISG15 shares 65.8% amino acid (aa) sequence identity with mouse ISG15. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ISG15 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.