

Anti-GAL4/LGALS4 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: PB9952-Fluoro594

About LGALS4

Galectin-4 is a protein that in humans is encoded by the LGALS4 gene. This gene is mapped to chromosome 19q13.2 based on an alignment of the LGALS4 sequence. The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. LGALS4 is an S-type lectin that is strongly underexpressed in colorectal cancer. The 323-amino acid LGALS4 protein contains 2 homologous, approximately 150-amino acid carbohydrate recognition domains and all amino acids typically conserved in galectins.

Overview

Product Name	Anti-GAL4/LGALS4 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P56470

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human GAL4, different from the related mouse and rat sequences by seven amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GAL4/LGALS4 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.