

## Anti-GAL4/LGALS4 Antibody Picoband® HRP Conjugated

Catalog Number: PB9952-HRP

### About LGALS4

Galectin-4 is a protein that in humans is encoded by the LGALS4 gene. This gene is mapped to chromosome 19q13.2 based on an alignment of the LGALS4 sequence. The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. LGALS4 is an S-type lectin that is strongly underexpressed in colorectal cancer. The 323-amino acid LGALS4 protein contains 2 homologous, approximately 150-amino acid carbohydrate recognition domains and all amino acids typically conserved in galectins.

### Overview

Product Name	Anti-GAL4/LGALS4 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P56470

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human GAL4, different from the related mouse and rat sequences by seven amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GAL4/LGALS4 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.