

Anti-MGP Antibody Picoband® Fluoro647 Conjugated

Catalog Number: PB9954-Fluoro647

About MGP

Matrix Gla protein (MGP) is an 84-residue vitamin K-dependent protein initially isolated from bovine bone. In addition, MGP is a 10-kD protein produced and secreted by vascular smooth muscle cells and chondrocytes and significantly accumulated in bone, cartilage, and dentin. It is also expressed at high levels in heart, kidney, and lung and is upregulated by vitamin D in bone cells. MGP has a high affinity binding to calcium ions, similar to other Gla-containing proteins. The protein acts as an inhibitor of vascular mineralization and plays a role in bone organization. And this gene is mapped to 12p12.3.

Overview

Product Name	Anti-MGP Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P08493

Technical Details

Immunogen	E.coli-derived human MGP recombinant protein (Position: Y20-F96). Human MGP shares 84.4% and 80.5% amino acid (aa) sequence identity with mouse and rat MGP, respectively.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MGP Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.