

Anti-Carboxypeptidase B2/CPB2 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: PB9997-Fluoro550

About CPB2

Carboxypeptidase B2 (CPB2), also known as carboxypeptidase U (CPU), plasma carboxypeptidase B (pCPB) or thrombin-activatable fibrinolysis inhibitor (TAFI), is an enzyme that, in humans, is encoded by the gene CPB2. CPB2 is synthesized by the liver and circulates in the plasma as a plasminogen-bound zymogen. When it is activated by proteolysis at residue Arg92 by the thrombin/thrombomodulin complex, CPB2 exhibits carboxypeptidase activity. Activated CPB2 reduces fibrinolysis by removing the fibrin C-terminal residues that are important for the binding and activation of plasminogen.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Carboxypeptidase B2/CPB2 Antibody Picoband® Fluoro550 Conjugated |
| Reactive Species | Human |
| Application | Recommended applications are based on the parent unconjugated antibody (IHC, WB, ELISA (Cap)). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q96IY4 |

Technical Details

| | |
|------------------|--|
| Immunogen | E. coli-derived human CPB2 recombinant protein (Position: K166-D388). Human CPB2 shares 86.1% and 85.7% amino acid (aa) sequence identity with mouse and rat CPB2, respectively. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm |

Suggested Dilutions

Optimal dilutions should be determined by end users.

1 Publications Citing This Product

1. PubMed ID: 19018804, INVOLVEMENT OF PROLYLCARBOXYPEPTIDASE IN THE EFFECT OF RUTAECARPINE ON THE REGRESSION OF MESENTERIC ARTERY HYPERTROPHY IN RENOVASCULAR HYPERTENSIVE RATS

Visit bosterbio.com/anti-cpb2-picoband-trade-antibody-pb9997-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Carboxypeptidase B2/CPB2 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.