

Human ANGPT2 Recombinant Protein

Catalog Number: PR00370-5

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human ANGPT2 Recombinant Protein
Description	This product is widely applicable in biological and medical research. It can serve as an antigen for antibody production, enabling the generation of specific antibodies following immunization.
Size	100 ug
Tag	
Form	Liquid
Source	
Formulation	Human ANGPT2 Recombinant Protein 0.7%, BSA 0.5%, glycerol 50%, Water 48.8%

Storage

Maintain refrigerated at 2–8°C for up to 2 weeks. For long-term storage, aliquot and store at -20°C to avoid repeated freeze-thaw cycles.

Purity

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human ANGPT2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.