

Human Myosin Id Recombinant Protein

Catalog Number: PR09126-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human Myosin Id Recombinant Protein
Description	This product is widely applicable in biological and medical research. It can serve as an antigen for antibody production, enabling the generation of specific antibodies following immunization.
Size	100ul
Tag	
Form	Liquid
Source	
Formulation	Human Myosin Id Recombinant Protein 20.7%, BSA 0.5%, glycerol 50%, Water 48.8%

Storage

Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Purity

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human Myosin Id Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.