

FITC conjugated Bovine Serum Albumin

Catalog Number: PR2G6FC

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	FITC conjugated Bovine Serum Albumin
Description	This product is widely used in biological and biomedical research. For example, it can serve as an antigen for antibody production, and immunization with this product enables the generation of specific antibodies.
Size	0.5 ml
Tag	
Form	Liquid
Source	
Formulation	BSA: 75%, FITC: 4.5%, Water: 20.5%

Concentration

1mg/ml

Storage

For long term storage, the antibody should be stored at -20°C and is stable for 12 months from the date of receipt. Avoid repeated freeze-thaw cycles. For short term storage, keep at 4°C for a week should be fine.

Purity

Protein A affinity chromatography

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[FITC conjugated Bovine Serum Albumin](#)

For Research Use Only. Not for use in diagnostic procedures.