

Human Vinculin Recombinant Protein

Catalog Number: PR9131

Overview

Product Name	Human Vinculin Recombinant Protein
Description	This product is widely applicable in biological and medical research. It can serve as an antigen for antibody production, enabling the generation of specific antibodies following immunization.
Size	0.2ml
Tag	
Form	Liquid
Source	
Formulation	Human Vinculin Recombinant Protein 20.7%, BSA 0.5%, glycerol 50%, Water 48.8%

Storage

Purity

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Vinculin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.