

Mouse IL-1 beta Recombinant Protein

Catalog Number: PRO401N

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse IL-1 beta Recombinant Protein
Description	This product is widely used in various biological and medical experiments. For example, it can be used in antibody preparation experiments as an antigen, and corresponding antibodies can be obtained after immunization.
Size	100ug
Tag	
Form	Lyophilized
Source	Polyclonal Goat IgG
Formulation	0.1% Mouse IL-1 beta Recombinant Protein, 5%Tris, 30% glycerol, 10% NaCl, 0.2% Dithiothreitol, 54.7% Water

Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
12 months from date of receipt, -20 to -70 °C as supplied.
1 month, 2 to 8 °C under sterile conditions after reconstitution.
6 months, -20 to -70 °C under sterile conditions after reconstitution.

Purity

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Mouse IL-1 beta Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.