

Human TNF alpha Recombinant Protein

Catalog Number: PROA21A

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human TNF alpha Recombinant Protein
Description	This product is widely used in various biological and medical experiments. For example, it can be used in antibody preparation experiments as an antigen, which can induce the production of corresponding antibodies after immunization.
Size	100ug
Tag	
Form	Lyophilized
Source	Polyclonal Goat IgG
Formulation	Human TNF alpha Recombinant Protein: 0.7%,BSA 0.5%, Trehalose: 17.7%, NACL: 30%, glycerol 51.1%

Storage

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Purity

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human TNF alpha Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.