

Recombinant Influenza A H7N9 Hemagglutinin/HA1 protein

Catalog Number: PROTA0A024E3P0

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Recombinant Influenza A H7N9 Hemagglutinin/HA1 protein
Description	Recombinant Influenza A H7N9 Hemagglutinin/HA1 protein expressed in Baculovirus with His-Tag. Sequence domain: 19-339aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Influenza A, A0A024E3P0, 19-339aa; DKICLGH HAVSNGTKVN TLTERGVEVV NATETVERTN IPRICKSGKR TVDLGQCGLL GTITGPPQCD QFLEFSADLI IERREGSDVC YPGKFVNEEA LRQILRESGG IDKEAMGFTY SGIRTNGATS ACRRSGSSFY AEMKWLLSNT DNAAFPQMTK SYKNTRKSPA LIVWGIHHSV STAEQTKLYG SGNKLVTGVS SNYQQSFVPS PGARPQVNGL SGRIDFHWLM LNPNDVTFTS FNGAFIAPDR ASFLRGKSMG IQSGVQVDAN CEGDCYHSGG TIISNLPFQN IDRAVGKCP RYVKQRSLLL ATGMKNVPEI PKGR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Recombinant Influenza A H7N9 Hemagglutinin/HA1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.