

## Mouse recombinant IL-17D (Interleukin-17D) protein, AF

Catalog Number: PROTA0A0B4J1G4

### Overview

Product Name	Mouse recombinant IL-17D (Interleukin-17D) protein, AF
Description	Interleukin 17D (IL-17D) belongs to the IL-17 family of cytokines, predicts a molecular mass of 20 kDa. IL-17 family is closely linked to host defense and immune response, and IL-17D is a novel cytokine in the IL-17 family of cytokines that has not been extensively investigated. It is highly secreted by fibrosarcoma tumor cells; in addition, ectopic expression of IL-17D in tumor cells recruits natural killer cells via the CCL2 production of endothelial cells.
Size	5ug,20ug,100ug
Tag	His Tag (N-term)
Form	Lyophilized
Source	Escherichia coli
Formulation	The protein was lyophilized from a 0.2 um filtered solution containing 20 mM sodium citrate, 0.2 M NaCl, pH 4.5. If you have any concerns or special requirements, please confirm with us.

### Storage

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

### Purity

>98% as determined by SDS-PAGE.

### Molecular weight

The protein has a calculated MW of 20.74 kDa. The protein migrates about 25 kDa under reducing condition (SDS-PAGE analysis).

### Activity

Testing in process

### Endotoxin

<0.1 EU per 1 ug of the protein by the LAL method.

### Amino Acid Sequence

ALRTGRRPARPRDCADRPEELLEQLYGRLAAGVLSAFHHTLQLGPREQARNASCPAGGRAADRRFRPPTNLRVSPWAYRISYDPARFPRYLPEAYCLCR  
GCLTGLYGEEDFRFRSTPVFSPAVVLRRTAACAGGRSVYAEHYITIPVGCTCVPEPDKSAZDSANSSMDKLLLGPADRPAGR with polyhistidine tag at  
the N-terminus.

### Reconstitution

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H<sub>2</sub>O to a concentration not less than 100 ug/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant IL-17D (Interleukin-17D) protein, AF

For Research Use Only. Not for use in diagnostic procedures.