

## CXorf49B (NM\_001145139) Human Recombinant Protein

Catalog Number: PROTA8MYA2

### Overview

|              |   |
|--------------|---|
| Product Name | CXorf49B (NM_001145139) Human Recombinant Protein                           |
| Description  | Recombinant protein of human chromosome X open reading frame 49B (CXorf49B) |
| Size         | 20 µg   |
| Tag          | C-Myc/DDK   |
| Form         | Frozen Solution in PBS Buffer   |
| Source       | HEK293T   |
| Formulation  | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol                        |

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MSSPDKVSVCGAGFDLEGGKKAGSRTASPGAPGAHSHGLDLGVPGSGDGKSESGFTDPEGFSFESESELIEQGRVVLWGREGRPGTPVDDQGDVVDYS  
FYLADEPAAIVPPPSVQGHFPPEGAAAEGSAENWADAEVGPGRDVLGHSPGKWQQASAGRLHLCGPGPVRAWKNPERGSKSRWSLRVDPQQPSAK  
GPTRLPETHSDSADSSDLPLMKVIGICRNEGSQAKPGSPKKRADTSRQASFHCKESYLPVGRFLTSAPRGLTPVAERP AVGELEDSPQKKMQSRAWG  
KVEVRPSCSGAAAAGALPQGLSRRKMAGGKSLGGASQLALGRGFPACGERLSAAPPEPATFPFSGVVRPQGMSSKPKPKHSSPGKKPAGRKTRESQ  
AAAREDNDPNRDEVPRACLPTHRPGLPRLSVRRGEFSSSDPNIRAPQLPGTSEPSAYSPGGLVPRRHAPSGNQPPVHPPRPERQQQPPGAQGCPRCI  
WLQREIEDLTQQLAAMQFLTDKFQDL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



CXorf49B (NM\_001145139) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.