

influenza a-h3n2 (Human) Hemagglutinin Recombinant Protein

Catalog Number: PROTC6KNH7

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	influenza a-h3n2 (Human) Hemagglutinin Recombinant Protein
Description	influenza a-h3n2 (Human) Hemagglutinin Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 17-345aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Influenza A-H3N2 (human), C6KNH7, 17-345aa; QKLPGN DNSTATLCLG HHAVPNGTIV KTITNDQIEV TNATELVQSS STGEICDSPH QILDGKNCTL IDALLGDPQC DGFQNKKDWL FVERSKAYSN CYPYDVPDYA SLRSLVASSG TLEFNNEFN WTGVTQNGTS SACIRRSKNS FFSRLNWLTH LNFKYPALNV TMPNNEQFDK LYIWGVHHPG TDKDQIFLYA QASGRITVST KRSQQTVPSPN IGSRPRVRNI PSRISIWWTI VKPGDILLIN STGNLIAPRG YFKIRSGKSS IMRSDAPIGK CNSECITPNG SIPNDKPFQN VNRITYGACP RYVKQNTLKL ATGMRNVPEK QTR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[influenza a-h3n2 \(Human\) Hemagglutinin Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.