

Human HSBP1L1 Recombinant Protein

Catalog Number: PROTC9JCN9

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human HSBP1L1 Recombinant Protein
Description	Human HSBP1L1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-74aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, C9JCN9, 1-74aa; MGSSHHHHHH SGLVPRGSH MGSM DVRGPE APGGRALRDA AENLFQELQE HFQALTATLN LRMEEMGNRI
EDLQKNVNDL MVQAGIENSI KEQMLKT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human HSBP1L1 Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.