

Human recombinant Galectin-9 protein, AF

Catalog Number: PROTO00182-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Galectin-9 protein, AF
Description	Galectin-9 (Gal-9) is a lectin family member and is one of the tandem repeat-type galectins containing two carbohydrate recognition domains (CRD) connected by a linker region in a single peptide chain. The CRD is responsible for beta-galactoside binding, and several binding partners for galectin-9 have been identified, including CD44, FGL, Gcgr, GLUT- 2, IgE, IgM, PDI, and Tim-3. Galectin-9 is constitutively presented in the small intestine, liver, uterine epithelial cells, skin epidermis, and esophageal epithelium. Galectin-9 contributing to cell growth, differentiation, adhesion, communication, and death. Accumulated evidence indicates that Gal-9 blockade promotes T cell immunity against tumors, suggesting that galectin-9 is a promising target for immunotherapy.
Size	5ug,20ug,100ug
Tag	His Tag (N-term)
Form	Lyophilized
Source	Escherichia coli
Formulation	The protein was lyophilized from a 0.2 um filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

Storage

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Purity

>98% as determined by SDS-PAGE.

Molecular weight

The protein has a calculated MW of 36.7 kDa. The protein migrates as 37 kDa under reducing condition (SDS-PAGE analysis).

Activity

Measured by its ability of the immobilized protein to support the adhesion of Jurkat cells. The ED₅₀ for this effect is <3 ug/mL.

Endotoxin

<0.1 EU per 1 ug of the protein by the LAL method.

Amino Acid Sequence

AFSGSQAPYLSPAVPFSGTIQGGLQDGLQITVNGTVLSSSGTRFAVNFQTGFSGNDIAFHFNPRFEDGGYVVCNTRQNGSWGPEERKTHMPFQKGMPF
DLCFLVQSSDFKVMVNGILFVQYFHRVPFHRVDTISVNGSVQLSYISFQPPGVWPANPAPITQTVIHTVQSAPGQMFSTPAIPPMYPHPAYPMPFITTI
LGGLYPSKSILLSGTVLPSAQRFHINLCSGNHIAFHLNPRFDENAVVRNTQIDNSWGSEERSLPRKMPFVRGQSFSVWILCEAHCLKVAVDGQHLEYYH
RLRNLPTINRLEVGGDIQLTHVQT with polyhistidine tag at the N-terminus.

Reconstitution

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H₂O to a concentration not less than 100 ug/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human recombinant Galectin-9 protein, AF

For Research Use Only. Not for use in diagnostic procedures.