

Human recombinant Podocalyxin protein

Catalog Number: PROTO00592

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Podocalyxin protein
Description	Human recombinant Podocalyxin protein expressed in HEK293 with His-Tag. Sequence domain: 23-429aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol.

Concentration

1 mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, O00592, 23-429aa; SPSPSPS PSQNATQTTT DSSNKTAPTP ASSVTIMATD TAQQSTVPTS KANEILASVK ATT LGVSSDS PGTTTTLAQQV
SGPVNTTVAR GGGSGNPSTT IESPKSTKSA DTTT VATSTA TAKPNTTSSQ NGAEDTTNSG GKSSHSVTTD LTSTKAEHLT TPHPTSPLSP
RQPTSTHPVA TPTSSGHDHL MKISSSSSTV AIPGYFTSP GMTTTLPSV ISQRTQQTSS QMPASSTAPS SQETVQPTSP ATALRTPTLP ETMSSSPTAA
STTHRYPKTP SPTVAHESNW AKCEDLETQT QSEKQLVLNL TGNTLCAGGA SDEKLISLIC RAVKATFNPA QDKCGIRLAS VPGSQTVVVK

EITHTKLPA KDVYERLKDK WDELKEAGVS DMKLGDQGPP EEAEDRFSMP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant Podocalyxin protein](#)

For Research Use Only. Not for use in diagnostic procedures.