

Rat recombinant UNC5H2/UNC5B protein

Catalog Number: PROTO08722

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Rat recombinant UNC5H2/UNC5B protein
Description	Rat recombinant UNC5H2/UNC5B protein expressed in HEK293 with hlgG-His-Tag. Sequence domain: 27-373aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Rat, O08722, 27-373aa; GIDSGG QALPDSFPSA PAEQLPHFLL EPEDAYIVKN KPVELHCRAF PATQIYFKCN GEWVSQKGHV TQESLDEATG
LRIREVQIEV SRQQVEELFG LEDYWCQCVA WSSSGTTKSR RAYIRIAYLR KNFDQEPLAK EVPLDHEVLL QCRPPEGVPV AEVEWLKNE
VIDPAQDTNF LLTIDHNLII RQARLSDTAN YTCVAKNIVA KRRSTTATVI VYVNGGWSSW AEWSPCSNRC GRGWQKRTRT CTNPAPLN
AFCEGQACQK TACTTVCPVD GAWTEWSKWS ACSTECAHWR SRECMAPPPQ NGGRDCSGLT LDSKNCTDGL CVLNQRTLND PKSRPLEPSG D

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Rat recombinant UNC5H2/UNC5B protein](#)

For Research Use Only. Not for use in diagnostic procedures.