

## Mouse Ogg1 Recombinant Protein

Catalog Number: PROTO08760

### Overview

Product Name	Mouse Ogg1 Recombinant Protein
Description	Mouse Ogg1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-345aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, O08760, 1-345aa; MGSSHHHHHH SGLVPRGSH MGSMLFRSWL PSSMRHRTLSSPALWASIP CPRSELRLDL VLASGQSFRW  
KEQSPAHWWSG VLADQVWTLT QTEDQLYCTV YRGDSSQVSR PTLEELET LH KYFQLDVSLA QLYSHWASVD SHFQRVAQKF QGVRLLRQDP  
TECLFSFICS SNNNIARITG MVERLCQAFG PRLIQLDDVT YHGFPNLHAL AGPEAETHLR KLGLGYRARY VRASAKAILE EQGGPAWLQQ  
LRVAPYEEAH KALCTLPGVG AKVADCICLM ALDKPQAVPV DVHVWQIAHR DYGWHPKTSQ AKGPSPLANK ELGNFFRN LW GPYAGWAQAV  
LFSADLRQPS LSREPPAKRK KGSKRPEG

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Ogg1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.