

## Human TRDMT1 Recombinant Protein

Catalog Number: PROTO14717

### Overview

Product Name	Human TRDMT1 Recombinant Protein
Description	Human TRDMT1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-391aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O14717, 1-391aa; MEPLRV LELYSVGGGM HHALRESCIP AQVVAIDVN TVANEVYKYN FPHTQLLAKT IEGITLLEFD RLSFDMILMS PPCQPFTRIG RQGDMTDSRT NSFLYILDIL PRLQKLPKYI LLENVKGFEV SSTRDLLIQT IENCGFYQE FLLSPTSLGI PNSRLRYFLI AKLQSEPLPF QAPGQVLMF PKIESVHPQK YAMDVENKIQ EKNVEPNISF DGSIQCSGKD AILFKLETAE EIHRKNQQDS DLSVKMLKDF LEDD TDVNQY LLPPKSLRLY ALLLDIVQPT CRRSVCFTKG YGSYIEGTGS VLQTAEDVQV ENIYKSLTNL SQEEQITKLL ILKLRVFTPK EIANLLGFPP EFGFPEKITV KQRYRLLGNS LNVHVAKLI KILYE

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human TRDMT1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.