

Human TACI/TNFRSF13B Recombinant Protein

Catalog Number: PROTO14836

Overview

Product Name	Human TACI/TNFRSF13B Recombinant Protein
Description	Human TACI/TNFRSF13B Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 1-165aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O14836, 1-165aa; ADLMSGLGRS RRGGRSRVDQ EERFPQGLWT GVAMRSCPEE QYWDPLLGTG MSCKTICNHQ SQRTCAAFCR
SLSCRKEQ GK FYDHLLRDCI SCASICGQHP KQCA YFCENK LRSPVNL PPE LRRQRSGEVE NNSD NSGRYQ GLEHRGSEAS PALPGLKLSA
DQVALVYSLE PKSCDKTHTC PPCPAPELLG GPSVFLFPPK PKDTLMISRT PEVTCVVVDV SHEDPEVKFN WYVDGVEVHN AKTKPREEQY
NSTYRVVSVL TVLHQDWLNG KEYKCKVSNK ALPAPIEKTI SKAKGQPREP QVYTLPPSRD ELTKNQVSLT CLVKG FYPSD IAVEWESNGQ
PENNYKTPP VLDS DGSFFL YSKLTVDKSR WQQGNV FSCS VMHEALHNHY TQKSLSLSPG KHHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human TACI/TNFRSF13B Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.