

Human Phospholipid Scramblase 1/PLSCR1 Recombinant Protein

Catalog Number: PROTO15162

Overview

Product Name	Human Phospholipid Scramblase 1/PLSCR1 Recombinant Protein
Description	Human Phospholipid Scramblase 1/PLSCR1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-288aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 30% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, O15162, 1-288aa; MDKQNSQMNA SHPETNLPVG YPPQYPPTAF QGPPGYSGYP GPQVSYPPPP AGHSGPGPAG FPVNPQPVYN QPVYNQPVGA AGVPWMPAPQ PPLNCPGLE YLSQIDQILI HQQIELLEVL TGFETNKNYE IKNSFGQRVY FAAEDTDCCT RNCCGSPRPF TLRIIDNMGQ EVITLERPLR CSSCCCPCCL QEIEIQAPPV VPIGYVIQTW HPCLPKFTIQ NEKREDVLKI SGPCVVCSCC GDVDFEIKSL DEQC VVGKIS KHWTGILREA FTDADNFGIQ FPLDL DVK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Phospholipid Scramblase 1/PLSCR1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.