

## Human Survivin Recombinant Protein

Catalog Number: PROTO15392

### Overview

Product Name	Human Survivin Recombinant Protein
Description	Human Survivin Recombinant Protein expressed in E. coli with Calmodulin Tag. Sequence domain: 1-142aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	Calmodulin Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O15392, 1-142aa; MGAPTLPP AWQPFLKDHR ISTFKNWPFL EGCCTPERM AEAGFIHCPT ENPDLAQCF FCFKELEGWE PDDDPIEEHK KHSSGCAFLS VKKQFEELTL GEFLKLDLDRER AKNKIAKETN NKKKEFEETA KKVRRRAIEQL AAMD

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Survivin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.