

## Human 8-oxoGuanine DNA Glycosylase/OGG1 Recombinant Protein

Catalog Number: PROTO15527

### Overview

Product Name	Human 8-oxoGuanine DNA Glycosylase/OGG1 Recombinant Protein
Description	Human 8-oxoGuanine DNA Glycosylase/OGG1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-345aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 40% glycerol

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O15527, 1-345aa; MPARALL PRRMGHRTLA STPALWASIP CPRSELRLDL VLPSGQSFRW REQSPAHWSG VLADQVWTLT QTEEQLHCTV YRGDKSQASR PTPDELEAVR KYFQLDVTLA QLYHHWGSVD SHFQEVAQKF QGVRLLRQDP IECLFSFICS SNNNIARITG MVERLCQAFG PRLIQLDDVT YHGFP SLQAL AGPEVEAHLR KLGLGYRARY VSASARAILE EQGGLAWLQQ LRESSYEEAH KALCILPGVG TKVADCICLM ALDKPQAVPV DVHMMWHIAQRDYSWHPTTSQ AKGSPQTNK ELGNFFRSLW GPYAGWAQAV LFSADLRQCR HAQEPPAKRR KGSKGPEG

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human 8-oxoGuanine DNA Glycosylase/OGG1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.