

Human FABP7/B-FABP Recombinant Protein

Catalog Number: PROTO15540

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human FABP7/B-FABP Recombinant Protein
Description	Human FABP7/B-FABP Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-132aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	25mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 2mM EDTA, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, O15540, 1-132aa; MVEAFCATWK LTNSQNFDEY MKALGVGFAT RQVGNVTKPT VIISQEGDKV VIRTLSFKN TEISFQLGEE
FDETTADDRN CKSVVSLDGD KLVHIQKWDG KETNFVREIK DGKMVMTLTF GDVVAVRHYE KA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human FABP7/B-FABP Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.