

TECK Thymus Expressed Chemokine Mouse Recombinant Protein (CCL25)

Catalog Number: PROTO35903

Overview

Product Name	TECK Thymus Expressed Chemokine Mouse Recombinant Protein (CCL25)
Description	TECK Mouse Recombinant produced in E. coli is a single, non-glycosylated, polypeptide chain containing 121 amino acids and having a molecular mass of 14.1kDa. The TECK is purified by proprietary chromatographic techniques.
Size	5ug, 20ug, 1mg
Tag	
Form	Sterile Filtered White lyophilized (freeze-dried) powder.
Source	Escherichia coli
Formulation	Filtered (0.2µm) and lyophilized from a concentrated (1mg/ml) solution in in 1×PBS, pH7.4.

Concentration

Storage

Lyophilized TECK although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution TECK should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Purity

Greater than 97.0% as determined by (a) Analysis by RP-HPLC and (b) Analysis by SDS-PAGE.

Amino Acid Sequence

QGAFEDCCLG YQHRIKWNVL RHARNYHQQE VSGSCNLRAV RFYFRQKVVC GNPEDMNVKR AIRILTARKR LVHWKSASDS QTERKKSNHM
KSKVENPNST SVRSATLGHP RMVMMPRKTN N

?

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



TECK Thymus Expressed Chemokine Mouse Recombinant Protein (CCL25)