

Human Kallikrein 10 Recombinant Protein

Catalog Number: PROTO43240

Overview

Product Name	Human Kallikrein 10 Recombinant Protein
Description	Human Kallikrein 10 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 34-276aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 1mM DTT, 30% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O43240, 34-276aa; ADPALLPQND TRLDPEAYGS PCARGSQPWQ VSLFNGLSFH CAGVLVDQSW VLTAAHCGNK PLWARVGDDH
LLLLQGEQLR RTRRSVVHPK YHQGSGPILP RRTDEHDLML LKLARPVVLG PRVRALQLPY RCAQPGDQCQ VAGWGTTAAR RVKYNKGLTC
SSITILSPKE CEVFYPGVVT NNMICAGLDR GQDPCQSDSG GPLVCDLTLQ GILSWGVPYPC GSAQHPAVYT QICKYMSWIN KVIRSNHHHH HH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Kallikrein 10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.