

Human Desert Hedgehog/Dhh (C23II) Recombinant Protein

Catalog Number: PROTO43323

Overview

Product Name	Human Desert Hedgehog/Dhh (C23II) Recombinant Protein
Description	Human Desert Hedgehog/Dhh (C23II) Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 23-198aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 0.15M NaCl, 10% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O43323, 23-198aa; MGSSHHHHHH SGLVPRGSH MGSMIIGPGR GPVGRRRYAR KQLVPLLYKQ FVPGVPERTL GASGPAEGRV
ARGSERFRDL VPNYNPDIIIF KDEENSGADR LMTERCKERV NALAIAMNM WPGVRLRVTE GWDEDGHHAQ DSLHYEGRAL DITTSRDRN
KYGLLARLAV EAGFDWVYYE SRNHVHVSVK ADNSLAVRAG G

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Desert Hedgehog/Dhh (C23II) Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.