

Human HtrA2/Omi (S306A) Recombinant Protein

Catalog Number: PROTO43464

Overview

Product Name	Human HtrA2/Omi (S306A) Recombinant Protein
Description	Human HtrA2/Omi (S306A) Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 134-458aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, O43464, 134-458aa; MAVSPPPAS PRSQYNFIAD VVEKTAPAVV YIEILDRHPF LGREVPISNG SGFVVAADGL IVTNAHVVD
RRRVRVRLLS GDTYEAVVTA VDPVADIATL RIQTKEPLPT LPLGRSADVR QGEFVVAMGS PFALQNTITS GIVSSAQRPA RDLGLPQTNV
EYIQTDAID FGNAGGPLVN LDGEVIGVNT MKVTAGISFA IPSDRLREFL HRGEKKNSSS GISGSQRRYI GVMMLTLSPS ILAELQLREP
SFPDVQHGVL IHKVLGSPA HRAGLRPGDV ILAIGEVMVQ NAEDVYEAVR TQSQLAVQIR RGRETLLYV TPEVTE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HtrA2/Omi (S306A) Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.