

Human KLF4 Recombinant Protein

Catalog Number: PROTO43474

Overview

Product Name	Human KLF4 Recombinant Protein
Description	Human KLF4 Recombinant Protein expressed in E. coli with His-Cam Tag. Sequence domain: 11-395aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Cam Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT, 0.1M NaCl.

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 80% by SDS-PAGE

Amino Acid Sequence

Human, O43474, 11-395aa; MAHHHHHHMA DQLTEEQIAE FKEAFSLFDK DGDGTITTKE LGTVMRSLGQ NPTEAELQDM INEVDADGNG TIDFPEFLTM MARKMKDTS EEEIREAFRV FDKDNGYIS AAELRHVMTN LGEKLTDEEV DEMIREADID GDGQVNYEEF VQMMTAKGSM AVSDALLPSF STFASGPAGR EKTLRQAGAP NNRWREELSH MKRLPPVLPG RPYDLAAATV ATDLESGGAG AACGGSNLAP LPRRETEEFN DLLDLDFILS NSLTHPPESV AATVSSSASA SSSSSPSSSG PASAPSTCSF TYPIRAGNDP GVAPGGTGGG LLYGRESAPP PTAPFNLADI NDVSPSGGFV AELLRPELDP VYIPPPQPQP PGGGLMGKFV LKASLSAPGS EYGSPSIVS SKGSPDGSHV VVVAPYNGGP PRTCPKIKQE AVSSCTHLGA GPPLSNGHRP AAHDFPLGRQ LPSRTTPTLG LEEVLSSRDC HPALPLPPGF HPHPGPNYPS FLPDQMOPV PPLHYQELMP PGSCMPEEPK PKRGRRSWPR KRTAT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human KLF4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.