

## Human PSTPIP1 Recombinant Protein

Catalog Number: PROTO43586

### Overview

Product Name	Human PSTPIP1 Recombinant Protein
Description	Human PSTPIP1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-416aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O43586, 1-416aa; MGSSHHHHHH SGLVPRGSH MGSMPQLQF KDAFWCRDFT AHTGYEVLLQ RLLDGRKMCK DMEELLRQRA  
QAEERYGKEL VQIARKAGGQ TEINSLRAS FDSLKQQMENV GSSHIQLALT LREELRSLEE FRERQKEQRK KYEAVMDRVQ KSKLSLYKKA  
MESKKTYEQK CRDADDAEQ A FERISANGHQ KQVEKSNKA RQCKDSATEA ERVYRQSI AQ LEKVRAEW EQ EHRTTCEAFQ LQEFDRILTIL  
RNALWVHSNQ LSMQCVKDDE LYEEVRLTLE GCSIDADIDS FIQAKSTGTE PPAPVPYQNY YDREVTPLTS SPGIQPSCGM IKRFSGLLHG  
SPKTTSLAAS AASTETLTPT PERNEG VYTA IAVQEIQGNP ASPAQEYRAL YDYTAQNPDE LDLSAGDILE VILEGEDGWW TVERNGQRGF  
VPGSYLEKL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human PSTPIP1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.