

ZNF207 (NM_001098507) Human Recombinant Protein

Catalog Number: PROTO43670

Overview

Product Name	ZNF207 (NM_001098507) Human Recombinant Protein
Description	Purified recombinant protein of Homo sapiens zinc finger protein 207 (ZNF207), transcript variant 3
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MGRK K K K Q L K P W C W Y C N R D F D D E K I L I Q H Q K A K H F K C H I C H K K L Y T G P G L A I H C M Q V H K E T I D A V P N A I P G R T D I E L I Y G M E G I P E K D M D E R R R L L E
Q K T Q E S Q K K K Q Q D D S D E Y D D D D S A A S T S F Q P Q P V Q P Q Q G Y I P P M A Q P G L P P V P G A P G M P P G I P L M P G V P P L M P G M P P V M P G M P P G L H H Q R K Y T
Q S F C G E N I M M P M G G M M P P G P G I P L M P G M P P G M P P V P R P G I P P M T Q A Q A V S A P G I L N R P P A P T A T V P A P Q P P V T K P L F P S A G Q M G T P V T S S T A S S
N S E S L S A S S K A L F P S T A Q A Q A A V Q G P V G T D F K P L N S T P A T T T E P P K P T F P A Y T Q S T A S T T S T N S T A A K P A A S I T S K P A T L T T T S A T S K L I H P D E D I S L E E R R A
Q L P K Y Q R N L P R P G Q A P I G N P P V G P I G G M M P P Q P G I P Q Q G M R P P M P P H G Q Y G G H H Q G M P G Y L P G A M P P Y Q G P P M V P P Y Q G G P P R P P M G M R P P
V M S Q G G R Y

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



ZNF207 (NM_001098507) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.