

## Human Siglec-6/CD327 Recombinant Protein

Catalog Number: PROTO43699

### Overview

Product Name	Human Siglec-6/CD327 Recombinant Protein
Description	Human Siglec-6/CD327 Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 27-347aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, O43699, 27-347aa; QERRFQL EGPESLTVQE GLCVLVPCL PTTLPASYYG YGYWFLEGAD VPVATNDPDE EVQEETRGRF HLLWDPRRKN  
CSLSIRDARR RDNAAYFFRL KSKWMKYGYT SSKLSVRVMA LTHRPNISIP GTLESGHPSN LTCSVPWVCE QGTPPIFSWM SAAPTSLGPR  
TTQSSVLTIT PRPQDHSTNL TCQVTFPGAG VTMERTIQLN VSYAPQKVAI SIFQGNAAF KILQNTSSLP VLEGQALRLL CDADGNPPAH  
LSWFQGFAL NATPISNTGV LELPQVGS AE EGDFTCRAQH PLGSLQISLS LFWHWKPEGR AGGV

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Siglec-6/CD327 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.