

Human SPOP Recombinant Protein

Catalog Number: PROTO43791

Overview

Product Name	Human SPOP Recombinant Protein
Description	Human SPOP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-374aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 50% glycerol, 0.2M NaCl, 2mM EDTA

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O43791, 1-374aa; MGSSHHHHHH SGLVPRGSH MSRVSPPPP AEMSSGPVAE SWCYTQIKVV KFSYMWTINN FSFCREEMGE VIKSSTFSSG ANDKLRWCLR VNPVKGLDEES KDYLSLYLLL VSCPKEVRA KFKFSILNAK GEETKAMESQ RAYRFVQGKD WGFKKFIRRD FLLDEANGLL PDDKLTLCFCE VSVVQDSVNI SGQNTMNMVK VPECRLADEL GGLWENSRT DCCLCVAGQE FQAHKAILAA RSPVFSAMFE HEMEEKKNR VEINDVEPEV FKEMMCFIYT GKAPNLDKMA DDLLAAADKY ALERLKMVCE DALCSNLSVE NAAEILILAD LHSADQLKTQ AVDFINYHAS DVLETSGWKS MVVSHPHLVA EAYRSLASQA CPFLGPPRKR LKQS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SPOP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.