

## Mouse recombinant CD200 protein

Catalog Number: PROTO54901-1

### Overview

Product Name	Mouse recombinant CD200 protein
Description	Mouse recombinant CD200 protein expressed in Baculovirus with His-Tag. Sequence domain: 31-232aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, O54901, 31-232aa; QVEVVTQDER KALHTTASLR CSLKTSQEPL IVTWQKKKAV SPENMVTYSK THGVVIQPAY KDRINVTELG  
LWNSSITFWN TTLEDEGCYM CLFNTEFGSQK VSGTACTLTY VQPIVHLHYN YFEDHLNITC SATARPAPAI SWKGTGTGIE NSTESHFHSN  
GTTSVTSILR VKDPKTQVGK EVICQVLYLG NVIDYKQSLD KG

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant CD200 protein

For Research Use Only. Not for use in diagnostic procedures.