

Mouse recombinant DKK1 protein

Catalog Number: PROTO54908

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant DKK1 protein
Description	Mouse recombinant DKK1 protein expressed in HEK293 with His-Tag. Sequence domain: 32-272aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	50mM MES buffer (pH 6.5) containing 30% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, O54908, 32-272aa; TLNSVLINSN AIKNLPPPLG GAGGQPGSAV SVAPGVLYEG GNKYQTLDN YQYPCAEDEE CGSDEYCSSP
SRGAAGVGGV QICLACRKRR KRCMRHAMCC PGNYCKNGIC MPSDHSFPR GEIESIEN LGNDHNAAAG DGYPRRTTLT SKIYHTKGQE
GSVCLRSSDC AAGLCCARHF WSKICKPVLK EGQVCTKHKR KGSHGLEIFQ RCYCGEGLAC RIQKDHQAS NSSRLHTCQR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant DKK1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.