

Mouse Nnmt Recombinant Protein

Catalog Number: PROTO55239

Overview

Product Name	Mouse Nnmt Recombinant Protein
Description	Mouse Nnmt Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-264aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, O55239, 1-264aa; MGSSHHHHHH SGLVPRGSH MGSHMESGFT SKDTYLSHFN PRDYLEKYYS FGSRHCAENE ILRHLLKNLF
KIFCLGAVKG ELLIDIGSGP TIYQLLSACE SFTEIIVSDY TDQNLWELQK WLKKEPGAFD WSPVVTVCD LEGNRMKGPE KEEKLRRAIK
QVLKCDVTQS QPLGGVSLPP ADCLLSTLCL DAACPDLPAY RTALRNLSL LKPGGFLVMV DALKSSYYMI GEQKFSSLPL GWETVRDAVE
EAGYTIEQFE VISQNYSTT SNNEGLFSLV GRKPGRSE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Nnmt Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.