

Human KIAA Recombinant Protein

Catalog Number: PROTO60268

Overview

Product Name	Human KIAA Recombinant Protein
Description	Human KIAA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-401aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, O60268, 1-401aa; METPEVP VGSLIDFGPE APTSSPLEAP PPVLQDGDGS LGDGASESET TESADSENDM GESPSHPSWD QDRRSSNES
FSSNQSTEST QDEETLALRD FMRGYVEKIF SGGEDLDQEE KAKFGEYCSSL ENGKGREWFA RYVSAQRCNS KCVSEATFYR LVQSFVAVLF
ECHQMDDFGP AKNLMTMCFT YYHIGKPQLL PPESREKPAQ SIDSYLKSN SWLAEKKDIA ERLKNTSAR TENVKGFFGG LETKLGKPLA
RRNEEDENKP QEKRPRAVTA YSPEDEKKG KIIYLYTHLKQ QPIWHTLRFW NAAFFDAVHC ERTKRSPTTR EKWCHMTQEE RDDSLRFNEN
ITFGQLGTFT HNMLAFGLNK KLCNDFLKKQ AVIGNLDEEQ YKLLSDHIEQ MATE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human KIAA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.