

Human SORBS3 Recombinant Protein

Catalog Number: PROTO60504

Overview

Product Name	Human SORBS3 Recombinant Protein
Description	Human SORBS3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-329aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, O60504, 1-329aa; MADGGSP FLGRRDFVYP SSTRDPSASN GGGSPARREE KKRKAARLKF DFQAQSPKEL TLQKGDIVYI HKEVDKNWLE GEHHGRLGIF PANYVEVLP ADEIPKPIKPP TYQVLEYGEA VAQYTFKGD L EVELSFRKGE HICLIRKVNE NWYEGRITGT GRQGIFPASY VQVSREPLR LCDDGQLPT SPRLTAAARS ARHPSSPSAL RSPADPIDLG GQTSRRTGF SFPTQEPRPQ TQNLGTPGPA LSHSRGPSHP LDLGTSSPNT SQIHWTPYRA MYQYRPQNE D ELELREGDRV DVMQQCDDGW FVGVSRR TQK FGTFPGNYVA PV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SORBS3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.