

Human NMT2 Recombinant Protein

Catalog Number: PROTO60551

Overview

Product Name	Human NMT2 Recombinant Protein
Description	Human NMT2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-498aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT, 0.1M NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, O60551, 1-498aa; MGSSHHHHHH SSGLVPRGSH MAEDSESAAS QQSLELDDQD TCGIDGDNEE ETEHAKGSPG GYLGAKKKKK KQKRKKEKPN SGGTKSDSAS DSQEIKIQQP SKNPSVPMQK LQDIQRAMEL LSACQGPARN IDEAAKHRYQ FWDTPVVKL DEVITSHGAI EPDKDNVRQE PYSLPQGMW DTLDSLDAEV LKELYLLNE NYVEDDDNMF RFDYSPEFLW WALRPPGWLL QWHCGVRVSS NKKLVGFISA IPANIRIYDS VKKMVEINFL CVHKLRSCR VAPVLIREIT RRVNLEGIFQ AVYTAGVVLV KPIATCRYWH RSLNPKLVE VKFSLSRNM TLQRTMKLYR LPDVTKTSGL RPMEPKDIKS VRELINTYLYL QFHLAPVMDE EEVAHWFLPR EHIIDTFVVE SPNGKLTDFL SFYTLPTVM HHPAHKSLKA AYSFYNIHTE TPLLDLMSDA LILAKSKGFD VFNALDLMEN KTFLEKLFKFG IGDGNLQYYL YNWRCPTGDS EKVGLVLQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human NMT2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.