

# Human Hematopoietic Prostaglandin D Synthase/HPGDS Recombinant Protein

Catalog Number: PROTO60760

## Overview

Product Name	Human Hematopoietic Prostaglandin D Synthase/HPGDS Recombinant Protein
Description	Human Hematopoietic Prostaglandin D Synthase/HPGDS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-199aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.15M NaCl, 1mM DTT

## Concentration

1mg/ml (determined by Bradford assay)

## Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

## Purity

> 95% by SDS-PAGE

## Amino Acid Sequence

Human, O60760, 1-199aa; MPNYKL TYFNMRGRAE IIRYIFAYLD IQYEDHRIEQ ADWPEIKSTL PFGKPILEV DGLTLHQSLA IARYLTKNTD LAGNTEMEQC HVDAIVDTLD DFMSCFWAE KKQDVKEQMF NELLTYNAPH LMQDLDTYLG GREWLIGNSV TWADFYWEIC STLLLVFKPD LLDNHPRLVT LRKKVQAIPA VANWIKRRPQ TKL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Hematopoietic Prostaglandin D Synthase/HPGDS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.